

CHEMISAFE®
SAFETY STORAGE CABINETS

MULTIRISK CSB SERIES SAFETY CABINETS FOR
NON-FLAMMABLE CHEMICALS



MULTIRISK CSB SERIES



SAFETY CABINETS FOR NON-FLAMMABLE CHEMICALS



CSB120



CSB120G



CSB60G



CSB120UB



CSB60



CSB60UB

MULTIRISK CSB SERIES

SAFETY CABINETS FOR NON-FLAMMABLE CHEMICALS



CSB124



CSB124G



CSB62



CSB62G



CSB63



CSB63G

MULTIRISK CSB SERIES



SAFETY CABINETS FOR NON-FLAMMABLE CHEMICALS

SPECIFICATIONS:

- Complying with the EN 61010 norm.
- Safety signs complying with ISO 3864 norm.
- Separation for "product category".
- Modular solutions.
- Completely manufactured from 1,5 mm thickness steel sheet, cold pressed, epoxy coated with acid resistant epoxy powders electro-statically applied at a temperature of 200 °C, which ensures a layer thickness higher than 90 microns on the whole surface of the cabinet.
- Basin form shelf with welded corners.
- Load capacity 80 kg for shelf, evenly distributed.
- Metal Hinged doors or with polycarbonate windows.
- Doors opening angle 110°.

CODE	OPTIONALS	EXTERNAL DIMENSIONS L x P x H mm
CSB60R	Powder coated shelf for ACIDS 60	590 x 500 x 30
CSB60RI	Stainless steel shelf for 60	590 x 500 x 30
CSB120R	Powder coated shelf for 120	1130 x 500 x 30
CSB120RI	Stainless steel shelf for 120	1130 x 500 x 30
TIM24	Daily Timer	
ALLARME	Alarm open	
CS125	Electro wired in metal box full with ON / OFF	230 x 230 x 250

CODE SOLID DOORS	CODE GLASSES DOORS	DESCRIPTION	EXTERNAL DIMENSIONS L x P x H mm	SHELVES - QUANTITY AND FINISH	CAPACITY BOTTLES 1 liter FOR CABINET	WEIGHT Kg
CSB60	CSB60G	Cabinet with 1 door - 1 compartment	600 x 600 x 1950	3 powder coated metal shelves	110 L	68
CSB62	CSB62G	Cabinet with 2 doors - 2 compartments	600 x 600 x 1950	ACIDS: 2 powder coated metal shelves BASES: 2 stainless steel shelves	110 L	70
CSB63	CSB63G	Cabinet with 3 doors - 3 compartments	600 x 600 x 1950	ACIDS: 2 powder coated metal shelves BASES: 2 stainless steel shelves TOXICS: 2 powder coated metal shelves	90 L	73
CSB120	CSB120G	Cabinet with 2 doors - 1 compartment	1140 x 600 x 1950	3 powder coated metal shelves	200 L	125
CSB124	CSB124G	Cabinet with 4 doors - 4 compartments	1140 x 600 x 1950	ACIDS: 2 powder coated metal shelves BASES: 2 stainless steel shelves TOXICS: 2 powder coated metal shelves HARMFUL: 2 powder coated metal shelves	200 L	135
CSB60UB	NO	Cabinet with 1 door - 1 compartment	600 x 600 x 700	2 powder coated metal shelves	40 L	35
CSB120UB	NO	Cabinet with 2 doors - 1 compartment	1140 x 600 x 700	2 powder coated metal shelves	80 L	50

Authorized dealer: