

**BOCE, BALON BOCE I KANISTRI
BOCA GRADUIRANA SA ŠIROKIM GRLOM**



KATALOŠKI BROJ	VOLUMEN ml	Ø mm	h mm	Ø GRILA mm
K00408	50	40	77	24
K00409	100	48	89	24
K00410	250	61	126	34,50
K00411	500	75	158	34,50
K00412	1000	93	200	45
K00413	2000	115	247	45

MATERIJAL:
BOCA: PE
ČEP: PP

BOCA GRADUIRANA ČETVRTASTA SA ŠIROKIM GRLOM



KATALOŠKI BROJ	VOLUMEN ml	DIMENZIJE mm	Ø GRILA mm
K00608	25	32 X 34 X 52	18
K00609	50	38 X 38 X 70	24
K00610	100	42 X 48 X 90	34,50
K00611	250	57 X 60 X 110	34,50
K00612	500	70 X 80 X 138	45
K00613	1000	80 X 103 X 176	58
K00614	2000	100 X 134 X 210	58

MATERIJAL:
BOCA: PE
ČEP: PP

BOCA GRADUIRANA USKO GRLO



KATALOŠKI BROJ	VOLUMEN ml	Ø mm	h mm	Ø GRILA mm
K00317	50	39	70	18,50
K00319	125	46	107	18,50
K00323	250	60	131	23
K00324	500	74	177	23
K00325	1000	94	205	33
K00326	2000	112	266	33

MATERIJAL:
BOCA: PE
ČEP: PP

BOCE GRADUIRANE ŠIROKO GRLO



KATALOŠKI BROJ	VOLUMEN ml	Ø mm	h mm	Ø GRILA mm	DIN GL
K01608	50	38	88	24	32
K01610	100	48	105	24	32
K01612	250	60	140	38	45
K01614	500	75	170	38	45
K01616	1000	95	206	55	63
K01618	2000	120	252	55	63

MATERIJAL:
BOCA I ČEP : PE

BOCE GRADUIRANE ŠIROKO GRLO



KATALOŠKI BROJ	VOLUMEN ml	Ø mm	h mm	Ø GRILA mm	DIN GL
K01621	50	38	88	24	32
K01623	100	48	105	24	32
K01625	250	60	140	38	45
K01627	500	75	170	38	45
K01629	1000	95	206	55	63
K01631	2000	120	252	55	63

MATERIJAL:
BOCA I ČEP : PP
MOGU SE AUTOKLAVIRATI

BOCA GRADUIRANA USKO GRLO



KATALOŠKI BROJ	VOLUMEN ml	Ø mm	h mm	Ø GRILA mm	DIN GL
K01595	50	38	92	13	18
K01597	100	48	108	13	18
K01599	250	60	150	19	25
K01601	500	75	182	19	25
K01603	1000	95	224	23	32

MATERIJAL:
BOCA I ČEP : PP
MOGU SE AUTOKLAVIRATI

BOCA USKO GRLO SA ČEPOM I SIGURNOSNIM PRSTENOM



KATALOŠKI BROJ	VOLUMEN ml	Ø mm	h mm	Ø GRILA mm
K01570	50	36	84	17,50
K01571	100	44	101	17,50
K01572	250	60	136	22,50
K01573	500	74	170	22,50
K01574	1000	90	222	22,50

MATERIJAL:
BOCA I ČEP : HDPE

BOCA ŠIROKO GRLO SA CRNIM ČEPOM I SIGURNOSNIM PRSTENOM



KATALOŠKI BROJ	VOLUMEN ml	Ø mm	h mm	Ø GRILA mm
K15710	250	60	137	30
K15711	500	74	177	30
K15712	1000	90	229	30

MATERIJAL:
BOCA I ČEP : HDPE

BOCA SA USKIM GRLOM SIVA



KATALOŠKI BROJ	VOLUMEN ml	Ø mm	h mm	Ø GRILA mm
K03171	50	39	70	18,50
K03191	125	46	107	18,50
K03231	250	60	131	23
K03241	500	75	158	23
K03251	1000	95	208	33
K03261	2000	112	288	33

MATERIJAL:
BOCA I ČEP : PE

BOCA ZA ISPIRANJE



KATALOŠKI BROJ	VOLUMEN ml	Ø mm	h mm
K00179	50	35	96
K00180	100	43	115
K00182	250	60	140
K00183	500	75	180
K00185	1000	95	220

MATERIJAL:
BOCA I ČEP : PE

BOCA ZA ISPIRANJE SA ŠIROKIM GRLOM GRADUIRANA



KATALOŠKI BROJ	VOLUMEN ml	Ø mm BOCE	Ø GRILA mm	DIN GL
K01637	250	60	38	45
K01638	500	75	38	45

MATERIJAL:
BOCA I ČEP : PE

BOCE SA ŠIROKIM GRLOM



KATALOŠKI BROJ	VOLUMEN ml	Ø mm BOCE	h mm	Ø GRILA mm
K01563	70	50	60	36
K01564	120	56	71	36
K01565	250	69	94	50
K01566	500	86	107	70
K01567	1000	111	128	85
K01568	2000	111	235	87
K01569	1500	111	182	87

MATERIJAL:
BOCA I ČEP :HDPE
CE
IVD

BOCA ŠIROKO GRLO GRADUIRANA

KATALOŠKI BROJ	VOLUMEN ml	DIMENZIJE mm	Ø GRILA mm	BOJA
K00479	2000	95X115X245	75	BIJELA
K00480	2000	95X115X245	75	SMEĐA

MATERIJAL:
BOCA I ČEP :HDPE
CE
IVD

BOCA BALON ŠIROKO GRLO

KATALOŠKI BROJ	VOL. Lit.	Ø GRILA UNUTARNJI mm	Ø GRILA VANJSKI mm	Ø BOCE mm	h mm
K01640	5	84	94	170	340
K01642	10	98	108	210	415

MATERIJAL:
BOCA I ČEP :HDPE

BOCA BALON USKO GRLO

KATALOŠKI BROJ	VOL. Lit.	Ø GRILA UNUTARNJI mm	Ø GRILA VANJSKI mm	Ø BOCE mm	h mm
K01644	5	52,50	62,50	170	340
K01646	10	52,50	62,50	210	425
K01648	25	79,50	95,50	280	565
K01650	50	79,50	95,50	350	700

MATERIJAL:
BOCA I ČEP :HDPE

BOCA BALON USKO GRLO SA PIPCEM



KATALOŠKI BROJ	VOL. Lit.	Ø GRILA UNUTARNJI mm	Ø GRILA VANJSKI mm	Ø BOCE mm	h mm
K01660	5	52,50	62,50	170	340
K01662	10	52,50	62,50	210	425
K01664	25	79,50	95,50	280	565
K01666	50	79,50	95,50	350	700

MATERIJAL:
BOCA I ČEP :HDPE

PIPAC 3/4"
KATALOŠKI BROJ: 375

KANISTAR SA PIPCEM



KATALOŠKI BROJ	VOLUMEN Lit.	DIMENZIJE mm
K01770	10	380x300x160
K01771	20	470x370x180

MATERIJAL:
KANISTER I ČEP :HDPE

INDUSTRIJSKI KANISTER



KATALOŠKI BROJ	VOLUMEN Lit.	DIMENZIJE mm	DIN GL
K01780	2	175x155x115	45
K01781	5	250x195x150	51
K01782	10	320x225x195	51
K01873	20	370x295x245	51

MATERIJAL:
KANISTER I ČEP :HDPE