

EPRUVETE (KEMIJSKE) SA RUBOM CILINDRIČNE

Borosilikatno staklo

Kataloški broj	Promjer (mm)	Duljina (mm)	Debljina stjenke(mm)	Pakiranje kom
ST 128799	10	60	1	100
ST 128800	10	100	1	100
ST 128900	13	100	1	100
ST 129000	13	130	1	100
ST 129100	16	100	1,2	100
ST 129200	16	160	1,2	100
ST 129300	18	180	1,2	100
ST 129310	18	200	1,2	100
ST 129320	20	200	1,2	100
ST 129330	20	200	1,2	100

