

## BOCE ZA KULTURU TKIVA 25 – 300 cm<sup>2</sup>



### STERILNE

DNA/RNA free

Dnase/RNase free

Pyrogens and endotoxins free

Cytotoxic substance free

KATALOŠKI BROJ	POVRŠINA RASTA cm <sup>2</sup>	VERZIJA	VOLUMEN ml	DIMENZIJE mm	MATERIJAL	PAKIRANJE/KOM
T 90025	25	VENT	3-8, max. 15	92 x 51 x 29	PS	36 X 10 / 360
T 90026	25	FILTER	3-8, max. 15	92 x 51 x 29	PS	36 X 10 / 360
T 90075	75	VENT	8-22, max. 65	155 x 87 x 40	PS	20 X 5 / 100
T 90076	75	FILTER	8-22, max. 65	155 x 87 x 40	PS	20 X 5 / 100
T 90150	150	VENT	15-45, max. 165	210 x 122 x 50	PS	12 X 3 / 36
T 90151	150	FILTER	15-45, max. 165	210 x 122 x 50	PS	12 X 3 / 36
T 90300	300	VENT	30-85, max. 410	275 x 170 x 50	PS	6 X 3 / 18
T 90301	300	FILTER	30-85, max. 410	275 x 170 x 50	PS	6 X 3 / 18

**BOCE ZA KULTURU TKIVA 115 – 150 cm<sup>2</sup>  
SA POKLOPCEM KOJI SE MOŽE VIŠE PUTA OTVARATI I ZATVARATI**



**STERILNE**

**DNA/RNA free**

**Dnase/RNase free**

**Pyrogens and endotoxins free**

**Cytotoxic substance free**

<b>KATALOŠKI BROJ</b>	<b>POVRŠINA RASTA cm<sup>2</sup></b>	<b>VERZIJA</b>	<b>VOLUMEN ml</b>	<b>DIMENZIJE mm</b>	<b>MATERIJAL</b>	<b>PAKIRANJE/ KOM</b>
T 90552	150	POKLOPAC	15-45, max. 100	210 X122X50	PS	6 X 3 / 18
T 90652	115	POKLOPAC / SA BARIJEROM	15-33, max. 100	210 X122X50	PS	6 X 3 / 18

## BOCE ZA KULTURU TKIVA 25 – 300 cm<sup>2</sup> SA PEEL-OF FOLIJOM



### STERILNE

DNA/RNA free

Dnase/RNase free

Pyrogens and endotoxins free

Cytotoxic substance free

KATALOŠKI BROJ	POVRŠINA RASTA cm <sup>2</sup>	VERZIJA	VOLUMEN ml	DIMENZIJE mm	MATERIJAL	PAKIRANJE/KOM
T 90028	25	PEEL-OF	3-8, max. 15	92 x 51 x 29	PS	36 X 10 / 360
T 90078	75	PEEL-OF	3-8, max. 15	92 x 51 x 29	PS	36 X 10 / 360
T 90153	150	PEEL-OF	8-22, max. 65	155 x 87 x 40	PS	20 X 5 / 100
T 90303	300	PEEL-OF	8-22, max. 65	155 x 87 x 40	PS	20 X 5 / 100
T 90653	115	PEEL-OF / SA BARIJEROM	15-45, max. 165	210 x 122 x 50	PS	12 X 3 / 36

**CLIPMAX 10 cm<sup>2</sup>**

KOMBINACIJA BOCE ZA KULTURU TKIVA I MIKROSKOPSKOG STAKLA 25 X 75 mm

**STERILNE**

DNA/RNA free

Dnase/RNase free

Pyrogens and endotoxins free

Cytotoxic substance free

KATALOŠKI BROJ	POVRŠINA RASTA cm <sup>2</sup>	VERZIJA	VOLUMEN ml	DIMENZIJE mm	MATERIJAL	PAKIRANJE/KOM
T 70010	10	FILTER, CLIP	-	93 x 23 x32	PS / COP	5 X 1 / 5



### ČEPOVI ZA BOCE ZA KULTURU TKIVA

Kataloški broj	Verzija	Dimenzija/mm	Materijal	Pak/kom Vrećica	Pak/kom Kutija
90825	VENT za T-25	24x16	PE	10	40
90826	FILTER za T-25	24x16	PE	10	40
90875	VENT za T-75	33x21	PE	10	40
90876	FILTER za T-75	33x21	PE	10	40
90850	VENT za T-150/300	39x24	PE	10	40
90856	FILTER za T_150/300	39x24	PE	10	40

### EPRUVETE ZA KULTURU TKIVA



**STERILNE**

**DNA/RNA free**


**Dnase/RNase free**

**Pyrogens and endotoxins free**

**Cytotoxic substance free**

<b>KATALOŠKI BROJ</b>	<b>POVRŠINA RASTA cm<sup>2</sup></b>	<b>VERZIJA</b>	<b>VOLUMEN ml</b>	<b>DIMENZIJE mm</b>	<b>MATERIJAL</b>	<b>PAKIRANJE/ KOM</b>
T 91106	20	VENT	5 - 8	16 X 122	PS	40 X 20 / 800
T 91243	10	KONIČNA- FILTER	3 - 10	30 X 105	PS	54 X 4 / 216

**STRUGAČ STANICA I STANIČNA ŠPATULA**




**STERILNE**

**DNA/RNA free**

**Dnase/RNase free**

**Pyrogens and endotoxins free**

**Cytotoxic substance free**

TIP	KATALOŠKI BROJ	DUŽINA mm	VERZIJA	ŠIRINA STRUGAĆA mm	MATERIJAL	PAKIRANJE / KOM
	T 99002	240	ROTIRAJUĆA	13	PE	150 X 1 / 150
	T 99003	300	ROTIRAJUĆA	20	PE	150 X 1 / 150
	T 99004	362	ROTIRAJUĆA	25	PE	100 X 1 / 100
	T 99010	195	FIKSNA	14	PE	150 X 1 / 100