

CRYO EPRUVETE, STALCI I KUTIJE




**STERILNE
SAMOSTOJEĆE
SA SILIKONSKIM RINGOM
ČEP NA NAVOJ VANJSKI
DNA/RNA free
Dnase/RNase free
Pyrogens and endotoxins free
Cytotoxic substance free**


KATALOŠKI BROJ	MAKSIMALNI VOLUMEN ml	PREPORUČENI VOLUMEN ml	DIMENZIJE mm	MATERIJAL	PAKIRANJE / KOM
T89012	0,9	0,7	17	PP	8 X 100
T89020	1,5	1,3	17	PP	8 X 100
T89040	3,5	3,3	17	PP	4 X 100
T89050	4,5	4,3	30	PP	4 X 100



STALAK ZA CRYO EPRUVETE



TIP	KATALOŠKI BROJ	TIP STALKA	KAPACITET / KOM	DIMENZIJE mm	MATERIJAL	PAKIRANJE / KOM
	T 99016	UNIVERZALNI, EPRUVETE OD 1,2 -5,00 ml	40	209 X 109 X 24	PP	1 / 10

KUTIJA ZA CRYO EPRUVETE



TIP	KATALOŠKI BROJ	CRYO EPRUVETE ml	KAPACITET / KOM	DIMENZIJE mm	MATERIJAL	PAKIRANJE / KOM
	T 99014	4,00 – 5,00	81	133 X 133 X 95	PP	1 / 10
	T 99015	1,20 - 2,00	81	133 X 133 X 45	PP	1 / 20