















TubeSpin® Bioreactors



TIP	KATALOŠKI BROJ	VERZIJA	VOLUMEN ml	DIMENZIJE mm	MATERIJAL	MAXIMUM RCF g	PAKIRANJE / KOM
	87015	15 ml konična	15	17,1 x 120	PP	15.500	20 X 40 / 800
	87017	15 ml cilindrična	15	17,1 x 120	PP	15.500	20 X 40 / 800
	87050	50 ml konična	50	30 x 115	PP	15.500	9 X 20 / 180
	787050	50 ml konična sa barijerom	50	30 x 115	PP	15.500	9 X 20 / 180
	87450	450 ml konična	450	96 x 145	PP	3.500	1 X 32 / 32
	87600	600 ml konična	600	100 x 183	PP	3.500	1 X 26 / 26
	99013	Stalak za 87600	3 x 87600	366x146 x70	PP		1 X 6 / 6
	99017	Stalak za 15 ml/ 50 ml	18 x 15 ml / 10 x 50 ml	209x98 x 62	PP		1 X 45 / 45
	87355	čep	450/600	70 x 21	PE		1 X 20 / 20
	87356	čep sa filterom	450/600	70 x 21	PE		1 X 20 / 20
	87414	Šejker ploča za 162/50 Kühner za 600 ml	162 x 50	800 x 420	ALU		1 X 1 / 1
	87611	Šejker ploča 32 Kühner za 600 ml	8 x 4	800 x 420	ALU		1 X 1 / 1

	87612	Šejker ploča 20 Kühner za 600 ml	5 x 4	500 x 420	ALU		1 X 1 / 1
	87613	Šejker ploča 16 Kühner za 600 ml	4 x 4	420 x 420	ALU		1 X 1 / 1
	87631	Šejker ploča 32 Infors Muttitron za 600 ml	8 x 4	850 x 470	ALU		1 X 1 / 1
	87633	Šejker ploča 16 Infors Minitron za 600 ml	4 x 4	480 x 420	ALU		1 X 1 / 1

TubeSpin® Bioreactors wit Septum



TIP	KATALOŠKI BROJ	VERZIJA	VOLUMEN ml	DIMENZIJE mm	MATERIJAL	MAXIMUM RCF g	PAKIRANJE / KOM
	86050	septum	50	30 x115	PP/PE	15.500	9 X 20/180