

## VAKUM FILTRACIJA "RAPID FILTERMAX"



### STERILNE

DNA/RNA free

Dnase/RNase free

Pyrogens and endotoxins free

Cytotoxic substance free

TIP	KATALOŠKI BROJ	POVRŠINA FILTERA cm <sup>2</sup>	VOLUMEN ml	TIP MEMBRANE μm	DIMENZIJE mm	MATERIJAL	PAKIRANJE / KOM
	T 99150	49	150	PES, 0,22	93X93X103	PS / PES	1 X 1 / 18
	T 99155	49	150	PES, 0,22	90X89X57	PS / PES	1 X 1 / 36
	T 99157		150		93X93X63	PS / PES	1 X 1 / 24
	T 99250	49	250	PES, 0,22	93X93X143	PS / PES	1 X 1 / 12
	T 99255	49	250	PES, 0,22	90X89X75	PS / PES	1 X 1 / 24
	T 99257		250		93X93X87	PS / PES	1 X 1 / 24
	T 99500	49	500	PES, 0,22	93X93X213	PS / PES	1 X 1 / 10
	T 99505	49	500	PES, 0,22	90X89X111	PS / PES	1 X 1 / 21
	T 99507		500		93X93X121	PS / PES	1 X 1 / 36
	T 99950	69	1000	PES, 0,22	111X111X285	PS / PES	1 X 1 / 9
	T 99955	69	1000	PES, 0,22	108X108X143	PS / PES	1 X 1 / 12
	T 99957		1000		111X111X160	PS / PES	1 X 1 / 20